

Project Background

Stroke is the 2nd leading causes of death, responsible for approximately 11% of total deaths respectively. *Up to 60% of stroke survivors suffer from upper extremity impairments.*



Stroke survivors continue to recover at home after inpatient rehabilitation ends. Progression may plateau without constant engagement in treatment.



■ Rehabilitation challenges faced by patients with hemiplegia

The restoration of distal motor function

The functional improvement of upper extremities primarily depends on the hand function. However, in the recovery process from upper extremity hemiplegia, the restoration of distal motor function comes later and is more strenuous than the restoration of proximal motor function.



Cognitive dysfunction

Cognitive impairment as evaluated with a comprehensive neuropsychological assessment is prevalent in stroke survivors even with successful clinical recovery. *Typically multiple domains and complex cognitive abilities are affected.*

A decreased mental state

Anxiety is a common psychological problem after stroke. Anxiety disorders or symptoms can also compromise rehabilitation; thus, the significance of patients' psychological status after stroke forms an essential element of their treatment process.

User Research and Interview



Huilan Zhou

Stroke Survivo

I feel:

Training is a heavy burden.

I want to:

Get rid of cumbersome training to return to normal life.



Changming Xu

Chief Physician, Renji Hospita

Patients feel:

Bored and painful in training.

Patients are:

Unwilling to carry out locomotor activities.

Q

Do patients after stroke have to undergo rehabilitation? Must be treated in hospital?

The post-stroke functional symptoms will not gradually recover voluntarily. It's impossible to stay in the hospital for a long time for rehabilitation after a stroke, so it's necessary to return home for rehabilitation.




Whether patients after stroke can fully recover to normal through rehabilitation training?

There is no way to fully restore the functional impairment after a stroke. Hard work in rehabilitation training can improve it, but it will get worse if patient don't work hard.



Hemiplegia is the paralysis of the muscles on one side of the body. The most common cause is stroke, which damages the corticospinal tracts in one hemisphere of the brain which extend from the lower spinal cord to the cerebral cortex.





Impaired hand function is one of the most frequently persisting consequences of stroke

When hemiplegia is caused by spinal cord injury , it may be called spinal hemiplegia.



The one side of the brain fails to control the muscles and other functions on the other side of the body.

Insight

Camera-based Systems



This method can be accurate but **involve** privacy issues. Also, the environment cannot be too cluttered, to avoid skeleton merging.

Tangible-interaction-objects-based System



These systems can reduce a learning times.


However, they are generally limited to a single training mode and have poor scalability.

Wearable-sensor-based System



EMG sensors have been widely utilized for rehabilitation systems. However, classification accuracy for stroke patients is typically much lower than that of healthy individuals.

We hypothesized that the proposed sensing and algorithm configuration could detect stroke patient hand gestures during the cognitive and motor tasks required in the serious games.



System Design

The SmartHand application will work with wearable devices to help users complete training. During training, the application will give action guidance.

At the same time, the camera module records hand movement data and transmits it back to the background in time. So that the algorithm can recommend more targeted daily activities and formulate rehabilitation plans next time.

Movement Record Intelligent Recommendation Rehabilitation training Information feedback module

Rehabilitation process feedback


Rehabilitation status update

Rehabilitation process visualization


Training feedback

(1=1) Training action correction

Data acquisition module



Current Rehabilitation



Storyboard

Daily time



Help users to complete hand movements through electrical stimulation.

Training time



Provide users interesting games, guide them to complete their training. Users would qualitatively be more enthusiastic about training.

After




Give users feedback about their recovery process. Regularly communicate with the doctor to update the rehabilitation plan.

■ Technical Research

Hardware


EMG detection and electrical stimulation are controlled through Arduino-based electrical stimulation module. The camera module is used to record gestures.





- Arduino UNO Board
- Electrical stimulation module
- Electrode slice
- @ Electrode slice
- 6 Camera
- 6 Power unit
- EMG sensor

Assembly and Functional Principles



EMG signals detection

The EMG sensor detects the myoelectric signal, and outputs EMS to promote the muscle contraction.




EMG signals detection


Apply electrical stimulation

Software and Coding

We use the Python-based MediaPipe to recognize gestures and compare them with standard gestures in the gesture library. It is converted into data by an algorithm, and the data is transmitted to the host.








Work Flow of Gesture Recognition

Mark hand coordinates

Template matching requires the establishment of a basic gesture library to compare the accuracy of patient hand movements.



Gesture Library

Cylinder grip (CG)


The FMA is an effective and detailed evaluation tool for assessing motor function after stroke. 11 upper extremity fine movements are selected from the FMA and are suitable for motor function rehabilitation. The movements include hand, wrist and forearm movements.



No motion (NM)








Opposition (O)

Hook like grasp (HG) Mass extension (ME)




Thumb adduction (TA)



Movements Generation Algorithm Framework

Spherical grip (SG)



Wrist volar flexion (WF)

1 Level of training

Assessed in the number of hand movements obtained by the vision module

² Preference of user Calculation of clicks from user training records and activity recommendation lists

Define the action requirements such as muscle training and action intensity corresponding to the training level

Prototyping and Test

Silicone Casting and Component Assembling





Assemblina


Welding





Front

Functional Model Test







Function test

Product Design

Camera Module



Ophthalmic lens Camera The camera module can be worn in a variety of ways, equipped with a gimbal to keep it stable.



Wearing performance

Wearable Module



A switch on the top of the device can be rotated to control the intensity of electrical stimulation.

The soft wrist strap wraps the arm to prevent slipping, and the magnetic buckle adjusts the tightness. The integrated electrodes containing hydrogel patches on the inside of the wristband conduct bidirectional current conduction.




Switch / Intensity knob EMG sensor Electrical stimulation module -Type-C port Silicone wristband - Air holes Electric conduction - EMS sensor

Interaction

How SmartHand helps people recover from hemiplegia?

SmartHand supports daily supplementary training through App. The algorithm builds a dynamic healing process for the user based on the data. App can give training recommendations and guidance, record and analyze data, and develop personalized training plans.



Architecture Diagram

The App assists patients in rehabilitation in the form of games. In addition to the game section, the App can view the rehabilitation progress and achievement wall to get the latest training plan and medical advice.



Interaction Prototype



Rehabitation tasks



User Profile

Serious Games

Two serious games were developed based on movement estimation. The games provide visual and audio feedback to the patients.

	Find the Zombie			Happy Salesman		
Group 1	MF	0	WD	TA	SG	FP
Group 2	ME	CG	WF	CG	FP	FS
Group 3	TA	SG	FP	HG	0	FS
Group 4	HG	NM	FS			

Patients select the correct target in the serious game and perform the corresponding movement. The game "Find the Zombie" was designed for both motor and cognitive function training. The game "Happy Salesman" was designed to train motor function and improve performance in ADLs.

Game 1: Find the Zombie



Complete the movement to flip the card

In this game, three cards appeared with a Zombie and two Plants on the front. Then, all cards were flipped over and randomly swapped positions. The patient was required to find which card is the one with the Zombie and perform the corresponding hand movement shown below that card.

Game 2: Happy Salesman

drumstick





Cover top

burger bread



burger bread





Coins

Customer's request Time left



Pass the right food to the customers



Level score



Daily use scene


The user owns a store that sells 7 types of food. Customers keep coming to the store to buy 1-3 types of food. Users need to pass the right food to the customers by performing the correct corresponding movement. The coins represents the number of objects that participants successfully "sold" to customers.



Game testing

Service System Mapping

SmartHand is equipped with a database developed by professional companies. User recovery data and professional medical advice can be uploaded to the database for users and hospitals. SmartHand has the function of continuous learning, calls the database in real time and dynamically adjusts rehabilitation plan.



Telemedicine to support recovery



Telemedicine can effectively solve the problem of shortage of medical resources. For long-term recovering patients, inperson medical care is often unnecessary Telemedicine improves efficiency and saves costs.

Data transfer and sharing between users



Build a user-centric database which data comes from doctors and users. The continuous learning of the algorithm is realized. Integrating usage information from others will improve the effect of rehabilitation.

■ Service Blueprint

